

## INSTRUCTION MANUAL FOR CHIRALPAK® IA-3, CHIRALPAK® IB-3, CHIRALPAK® IC-3, CHIRALPAK® ID-3, CHIRALPAK® IE-3 and CHIRALPAK® IF-3


### <Normal Phase>

**Please read this instruction sheet completely before using these columns**

**These columns can also be used in reversed phase mode.  
Please refer to the corresponding instruction sheet for details.**


#### Column Description

#### AMYLOSE-BASED



CHIRALPAK® IA-3	CHIRALPAK® ID-3	CHIRALPAK® IE-3	CHIRALPAK® IF-3
Amylose tris(3,5-dimethylphenylcarbamate) <b>immobilized on 3µm silica-gel</b>	Amylose tris(3-chlorophenylcarbamate) <b>immobilized on 3µm silica gel</b>	Amylose tris(3,5-dichlorophenylcarbamate) <b>immobilized on 3µm silica gel</b>	Amylose tris(3-chloro-4-methylphenylcarbamate) <b>immobilized on 3µm silica gel.</b>
R =	R =	R =	R =

#### CELLULOSE-BASED



CHIRALPAK® IB-3	CHIRALPAK® IC-3
Cellulose tris(3,5-dimethylphenylcarbamate) <b>immobilized on 3µm silica-gel</b>	Cellulose tris(3,5-dichlorophenylcarbamate) <b>immobilized on 3µm silica-gel</b>
R =	R =

Shipping solvent: **n-Hexane / alcohol solvent mixture**

All columns have been pre-tested before packaging. Test parameters and results, as well as the Column Lot Number, are included on a separate (enclosed) page.

## Operating Instructions

	<b>250 x 2.1 mm i.d.</b> <b>150 x 2.1 mm i.d.</b> Analytical columns	<b>50 x 4.6 mm i.d.</b> Analytical columns	<b>100 x 4.6 mm i.d.</b> Analytical columns	<b>150 x 4.6 mm i.d.</b> <b>250 x 4.6 mm i.d.</b> Analytical columns
Flow rate direction	As indicated on the column label			
Typical Flow rate	0.1 to 0.5 ml/min	0.5 to 5 ml/min	0.5 to 4 ml/min	0.5 to 2.5 ml/min
Temperature	0 to 40°C			

NOTES: The column is stable to HPLC pressures. At a given temperature, the column back pressure is linearly proportional to the flow rate.

## Method Development / Normal Phase

### A - Mobile phases

CHIRALPAK® IA-3, IB-3, IC-3, ID-3 and IE-3 can be used *with all ranges of organic miscible solvents*, progressing from the traditional mobile phases used with other DAICEL columns (mixtures of alkanes/alcohol, pure alcohol or acetonitrile (CH<sub>3</sub>CN)) to mobile phases containing methyl *tert*-butyl ether (MtBE), tetrahydrofuran (THF), dichloromethane (DCM), chloroform (CHCl<sub>3</sub>), ethyl acetate (EtOAc) among others.

### B - Method Development - Screening

When developing methods we would recommend a screening approach.

- The conditions described in Table 1 should be used as a **Primary Screening**.
- If the compound or compound series are not soluble in any of these mobile phases, we recommend progressing directly to the **Secondary Screening** (Table 2).

**Table 1. Immobilized Primary Screening Solvents**

<b>Primary solvent mixtures</b>	Alkane <sup>①</sup> /2-PrOH	Alkane <sup>①</sup> /EtOH	Alkane <sup>①</sup> /MtBE/EtOH <sup>②</sup>	Alkane <sup>①</sup> /THF <sup>③</sup>	Alkane/DCM <sup>④</sup> /EtOH
<b>Typical starting conditions</b>	80:20	80:20	0:98:2	70:30	50:50:2
<b>Advised optimisation range</b>	99:1 to 50:50	99:1 to 50:50	80:20:0 to 0:40:60	95:5 to 0:100	85:15:0 to 0:80:20

① Alkane = n-Hexane, iso-Hexane or n-Heptane. Some small selectivity differences may sometimes be found.

② In absence of alkane, methanol is more efficient than ethanol when combined with MtBE.

③ In the case of no environmental restrictions, **use of DCM is preferred to THF** in terms of better enantioselectivity that the former may induce.

④ For excessively retained samples, addition of ethanol up to 20% in pure DCM would be helpful.

If a suitable chiral separation is not found using the Immobilized Primary Screening strategy, we recommend an Immobilized Secondary Screening to be applied using the following conditions:

**Table 2. Immobilized Secondary Screening Solvents**

Secondary solvent mixtures	EtOAc <sup>①</sup> /Alkane <sup>②</sup>	CH <sub>3</sub> CN <sup>③</sup> /Alcohol <sup>④</sup>
Typical starting conditions	50:50	100:0
Advised optimisation range	20:80 to 100:0	100:0 to 0:100

- ① Alcohols (④) or THF can be added into EtOAc to enhance the eluting strength for strongly retained compounds.
- ② Alkane: n-Hexane, iso-Hexane or n-Heptane. Some small selectivity differences may sometimes be found.
- ③ Transfers between alkane mixtures and CH<sub>3</sub>CN are preferably made with a transition in alcohol in order to avoid miscibility issues.
- ④ Alcohol: MeOH, EtOH and 2-PrOH.

**Note:** All solvent proportions indicated in this manual are by volume.

### C – General Comments

- ⇒ Additional solvent combinations such as CHCl<sub>3</sub>/Alkane, 1,4-Dioxane/Alkane, Toluene/Alkane or Acetone/Alkane can also be investigated with CHIRALPAK® IA-3, IB-3, IC-3, ID-3, IE-3 and IF-3 columns.
- ⇒ The typical starting conditions represent the mobile phases of upper middle eluting strength. Under such conditions, most of the analytes can be eluted within a reasonable time range with a good probability of full resolution of the enantiomers.
- ⇒ Toluene, MtBE and chlorinated solvents can be used in their pure form as the mobile phase.
- ⇒ For fast eluting solvents, such as THF, we recommend to add alkane in order to modulate the retention.
- ⇒ Detection with a regular UV detector may become difficult depending on a combination of sample and mobile phase (e.g. EtOAc, high percentages of DCM). In these cases an alternative detector, such as RI detector or ELSD (Evaporative Light Scattering Detector), may be more effective than the UV detector.

### D – Additives

For basic or acidic samples, it is necessary to incorporate an additive into the mobile phase in order to optimise the chiral separation.

- ① It has been found that certain amines, such as EDA and AE induce much better behaviour for certain basic compounds than the most commonly used DEA.

☞ The addition of a low percentage of an alcohol (e.g. 2% EtOH or MeOH) in the mobile phase may be helpful to ensure the miscibility of EDA and AE with the low polarity mobile phases.

Basic Samples require Basic additives	Acidic Samples require Acidic additives
Diethylamine (DEA) 2-Aminoethanol (AE) <sup>①</sup> Ethylenediamine (EDA) <sup>①</sup> Butyl amine (BA)	Trifluoroacetic acid (TFA) Acetic acid Formic acid
< 0.5% Typically 0.1%	< 0.5% Typically 0.1%

⇒ **STRONGLY BASIC** solvent additives or sample solutions **MUST BE AVOIDED**, because they are likely to damage the silica gel used in this column.

## Column Care / Maintenance

- ❑ The use of a guard cartridge is highly recommended for maximum column life.
- ❑ Samples should be dissolved in the mobile phase. The mobile phase and the sample solution should be filtered through a membrane filter of approximately 0.5µm porosity to ensure that there is no precipitate before using.

### ☞ Column cleaning and regeneration procedures

Following extensive use of the column in multiple solvents there may be a change in column reproducibility. In order to ensure consistent performance, a regeneration method may be implemented to eliminate any change in chiral recognition due to the history of the column (mobile phases, additives...).

- Flush with ethanol at 0.5 ml/min for 30 min, followed by 100% THF at 0.5 ml/min for 2 hours.
- Flush with ethanol at **0.05 ml/min<sup>(\*)</sup>** for 300 min.

(\*) This low flow rate would be critical for the column performance.

If this is not successful, then try with 100% N,N-dimethylformamide (DMF) or N,N-dimethylacetamide (DMAC) at 0.3 ml/min for 3 hours instead of the THF flush.

⇒ This procedure is also recommended for switching between reversed phase and normal phase.

## Column storage

- For column storage, remove the acidic or basic additives by flushing the column with the same mobile phase without the additive. Columns can be stored with the additive-free mobile phases.

**Operating these columns in accordance with the guidelines outlined here will result in a long column life.**

⇒ If you have any questions about the use of these columns, or encounter a problem, please email [questions@chiraltech.com](mailto:questions@chiraltech.com) or call 800-6-CHIRAL for assistance.

Part Number	Name	Particle Size	Internal Diameter	Column Length	Product Type
80511	CHIRALPAK IA-3	3	4.0	10	Guard Cartridges (3)
80522	CHIRALPAK IA-3	3	4.6	50	Analytical
80523	CHIRALPAK IA-3	3	4.6	100	Analytical
80524	CHIRALPAK IA-3	3	4.6	150	Analytical
80525	CHIRALPAK IA-3	3	4.6	250	Analytical
80592	CHIRALPAK IA-3	3	2.1	50	Analytical
80593	CHIRALPAK IA-3	3	2.1	100	Analytical
80594	CHIRALPAK IA-3	3	2.1	150	Analytical
80595	CHIRALPAK IA-3	3	2.1	250	Analytical

Part Number	Name	Particle Size	Internal Diameter	Column Length	Product Type
81511	CHIRALPAK IB-3	3	4.0	10	Guard Cartridges (3)
81522	CHIRALPAK IB-3	3	4.6	50	Analytical
81523	CHIRALPAK IB-3	3	4.6	100	Analytical
81524	CHIRALPAK IB-3	3	4.6	150	Analytical
81525	CHIRALPAK IB-3	3	4.6	250	Analytical
81592	CHIRALPAK IB-3	3	2.1	50	Analytical
81593	CHIRALPAK IB-3	3	2.1	100	Analytical
81594	CHIRALPAK IB-3	3	2.1	150	Analytical
81595	CHIRALPAK IB-3	3	2.1	250	Analytical

Part Number	Name	Particle Size	Internal Diameter	Column Length	Product Type
83511	CHIRALPAK IC-3	3	4.0	10	Guard Cartridges (3)
83522	CHIRALPAK IC-3	3	4.6	50	Analytical
83523	CHIRALPAK IC-3	3	4.6	100	Analytical
83524	CHIRALPAK IC-3	3	4.6	150	Analytical
83525	CHIRALPAK IC-3	3	4.6	250	Analytical
83592	CHIRALPAK IC-3	3	2.1	50	Analytical
83593	CHIRALPAK IC-3	3	2.1	100	Analytical
83594	CHIRALPAK IC-3	3	2.1	150	Analytical
83595	CHIRALPAK IC-3	3	2.1	250	Analytical



Part Number	Name	Particle Size	Internal Diameter	Column Length	Product Type
84511	CHIRALPAK ID-3	3	4.0	10	Guard Cartridges (3)
84522	CHIRALPAK ID-3	3	4.6	50	Analytical
84523	CHIRALPAK ID-3	3	4.6	100	Analytical
84524	CHIRALPAK ID-3	3	4.6	150	Analytical
84525	CHIRALPAK ID-3	3	4.6	250	Analytical
84592	CHIRALPAK ID-3	3	2.1	50	Analytical
84593	CHIRALPAK ID-3	3	2.1	100	Analytical
84594	CHIRALPAK ID-3	3	2.1	150	Analytical
84595	CHIRALPAK ID-3	3	2.1	250	Analytical

Part Number	Name	Particle Size	Internal Diameter	Column Length	Product Type
85511	CHIRALPAK IE-3	3	4.0	10	Guard Cartridges (3)
85522	CHIRALPAK IE-3	3	4.6	50	Analytical
85523	CHIRALPAK IE-3	3	4.6	100	Analytical
85524	CHIRALPAK IE-3	3	4.6	150	Analytical
85525	CHIRALPAK IE-3	3	4.6	250	Analytical
85592	CHIRALPAK IE-3	3	2.1	50	Analytical
85593	CHIRALPAK IE-3	3	2.1	100	Analytical
85594	CHIRALPAK IE-3	3	2.1	150	Analytical
85595	CHIRALPAK IE-3	3	2.1	250	Analytical

Part Number	Name	Particle Size	Internal Diameter	Column Length	Product Type
86511	CHIRALPAK IF-3	3	4.0	10	Guard Cartridges (3)
86522	CHIRALPAK IF-3	3	4.6	50	Analytical
86523	CHIRALPAK IF-3	3	4.6	100	Analytical
86524	CHIRALPAK IF-3	3	4.6	150	Analytical
86525	CHIRALPAK IF-3	3	4.6	250	Analytical
86592	CHIRALPAK IF-3	3	2.1	50	Analytical
86593	CHIRALPAK IF-3	3	2.1	100	Analytical
86594	CHIRALPAK IF-3	3	2.1	150	Analytical
86595	CHIRALPAK IF-3	3	2.1	250	Analytical



**CHIRALPAK® IA-3**  
Analytical HPLC applications

CHIRALPAK® IA-3 (4.6x50mm)





Retention time (seconds)

Hexane-EtOH 80:20 (+0.1% AE)  
Flow rate: **5.0** ml/min, Temp.: 25°C



Retention time (seconds)

DCM-MeOH 98:2 (+0.1% AE)  
Flow rate: **5.0** ml/min, Temp.: 25°C





Retention time (seconds)

MeOH 100%  
Flow rate: **5.0** ml/min, Temp.: 25°C



**CHIRALPAK® IC-3**  
Analytical HPLC applications

CHIRALPAK® IC-3 (4.6x50mm)





Retention time (seconds)

MtBE-EtOH 98:2  
Flow rate: **5.0** ml/min, Temp.: 25°C



Retention time (seconds)

Hexane-EtOH 80:20 (+0.1% TFA)  
Flow rate: **5.0** ml/min, Temp.: 25°C



Retention time (seconds)

Hexane-IPA 90:10  
Flow rate: **6.0** ml/min, Temp.: 25°C