Chiral Technologies

Protecting Your Chiral Columns for Optimal Performance

Dr. Weston Umstead

Technology and Business Development Manager

Overview

- Pressure Issues
- Peak Shape Issues
- Retention or Selectivity Issues
- How can we help?



Issues with Pressure

ებეიტეი

აფითლიტი

ĊH₃

What Causes Pressure Issues?



- Blockage within the system
- Plugged inlet frit
 - Insoluble material in the sample
 - Particulate from wearing pump seals or rotors
 - Mobile phase incompatibility with sample solubility
 - Mobile phase component crashes out inside the column
 - Bacterial growth
- Mobile phase column dimension combinations



©2022 Daicel Chiral Technologies. All rights reserved.

Peak Shape Problems, Dolan, LCGC North America, 2008, pg. 610-616

Retention/Selectivity

(a)



What Causes Pressure Issues?

Figure 1: Peak splitting due to partial column blockage. Chromatograms: (a) sample and (b) analytical standard.

Figure 2: Characteristic peak shapes for (a, b) a partially blocked frit .

(b)



Figure 3: Schematic of flow path through column, where arrows represent sample before separation: (a) normal flow, (b) flow disturbance due to partially blocked frit,





How Can I Diagnose This Issues?



Reverse the column direction and inject a sample

- If the chromatography is restored, it's likely a blocked inlet frit
- If not, it's likely something else*
- QC the column and compare its performance to the test chromatogram that was provided with the column
- Try a different (larger i.d. or shorter length) column dimension.



Peak Shape Problems, Dolan, LCGC North America, 2008, pg. 610-616



Preemptive Pressure Issue Solutions

- Insoluble material in the sample
 - Filter sample with a 0.5 μm or 0.2 μm filter
- Particulate from wearing pump seals or rotors
 - Regular system maintenance
 - Use an in-line filter to catch particulate
- Mobile phase incompatibility with sample solubility
 - Ensure sample solvent/mobile phase match
 - Small solubility study to check the sample doesn't crash out







- Reversed Phase Short Term and Long Term
 - Okay to leave the column connected to the system (short term)
 - Flush the column with water/ACN (60:40) WITHOUT salts or buffer





Normal Phase or Polar Organic Short Term (less than a few days)

- Okay to leave the column connected to the system
- Flush the column with Ethanol
- Normal Phase or Polar Organic Long Term (week or longer)
 - Flush the column with Ethanol
 - Store in 90:10 Hex/IPA
 - If the column was used for Polar Organic mode (especially 100% ACN or 100% MeOH, ethanol flush is essential!)











- Protein Columns Short Term (less than a few days)
 - Flush the column with mobile phase that doesn't contain any salt or buffer (water/IPA 90:10)
 - Can be stored at ambient temperature in water/IPA
 - 85:15 for AGP and CBH
 - 90:10 for HSA/BSA/DSA/MSA/RSA
- Protein Columns Long Term (week or longer)
 - Same conditions as short term storage, except the column(s) should be stored in the refrigerator



After the Pressure Rises



- Particulate from wearing pump seals or rotors
 - Must replace the inlet frit
- Insoluble material from the sample or sample has crashed out of solution
 - Reverse flow of the column disconnect the detector
 - Find what the insoluble material is soluble in
- Most of the time a plugged inlet frit requires replacement of the column



Peak Shape Issues



800080 0

თეთთევ

ĊH₃

What Causes Peak Shape Issues?



- Plugged Inlet Frit
 - Poor sample cleanup or particulate from the instrument
- Voids in the packing
 - Mechanical or thermal "shock"
- Chemical modifications
 - Biological samples strongly adsorbing analytes or sample components
 - Column history or memory effect
- Won't discuss method development issues additives, sample overload, extra-column effects



Injection Solvent/Mobile Phase Mismatch





Practical HPLC Method Development, Snyder and Kirkland, 1997, pg. 221



©2022 Daicel Chiral Technologies. All rights reserved.

Peak Shape Problems, Dolan, LCGC North America, 2008, pg. 610-616







Figure 2: Characteristic peak shapes for (a, b) a partially blocked frit a

(c) void or column collapse.



(c) flow disturbance due to a column void.



(c)

What Causes Peak Shape Issues?



How Can I Diagnose This Issue?



- Plugged Inlet Frit
 - QC the column and compare to the test chromatogram
- Void in the packing
 - Open the column and remove the frit
- Chemical modifications
 - Try a different column or a new column
 - QC the column and compare to the test chromatogram



Preemptive Peak Shape Issue Solutions



- Plugged Inlet Frit
 - Filtering the sample good sample cleanup
 - Inline filter
 - Guard column
- Voids in the packing
 - Do not drop or hit the column on a hard surface
 - Operating at reasonable pressure and flow rate
 - Prep columns bent inlet frits
 - Ensure proper equilibration when changing conditions
 - Quick rotation of sample-injection valve



Preemptive Peak Shape Issue Solutions













Preemptive Peak Shape Issue Solutions







Plugged Inlet Frit Filtering the sample

Preemptive Peak Shape Issue Solutions

- Inline filter
- Guard column
- Voids in the packing
 - Do not drop or hit the column on a hard surface
 - Operating at reasonable pressure and flow rate
 - Prep columns bent inlet frits
 - Ensure proper equilibration when changing conditions
 - Quick rotation of sample-injection valve

Chemical modifications

- Good sample clean up
- Dedicate columns to specific conditions i.e. acidic OD-H, neutral OD-H, and a basic OD-H
- Dedicate columns to specific mobile phases i.e. once a coated phase is put in reversed phase mobile phase, it must stay reversed phase
 - Even immobilized phases shouldn't be moved back and forth on a regular basis







After Poor Peak Shape Occurs



- Plugged Inlet Frit
 - Washing the column in reverse
 - Replace the inlet frit
- Voids in the packing
 - Repack the column with new CSP
- Chemical modifications
 - Recondition or repack
 - Replace the column



Regeneration Examples







Regeneration Examples





Figure 1: Initial Performance Check of CHIRALPAK[®] IA-3 before Backflushing and Regeneration



Figure 2: Performance Check of CHIRALPAK[®] IA-3 after Backflushing and Regeneration





<section-header><section-header>

ĊНз



What Causes Changes in Retention or Selectivity?

- Adsorbed sample
 - Biological samples, petroleum-based samples
- Chemical Attack
 - Incompatible additives or solvents
 - Too high/low pH
- Won't discuss method development issues scaling from smaller to larger particles/columns, instrument differences, method robustness



How Can I Diagnose This Issue?



- Adsorbed sample
 - QC the column and compare to the test chromatogram
 - Open the column and inspect the CSP
- Chemical Attack
 - Presence of a white powdery solid exiting the column



Preemptive Retention or Selectivity Issue Solutions

R: N H

• Adsorbed sample

- Good sample clean up
 - SPE cartridges, achiral C18 column, membrane filtration, liquid extraction

Chemical Attack

- Develop methods with a pH high than 2, but lower than 8 (or 9 with the correct buffer)
 - Ligand cleavage or silica degradation
- Avoid harsh acidic additives like HCl or harsh basic additives like NaOH
 - Preferred options for acids = TFA, Formic acid, Acetic acid
 - Preferred options for bases = DEA, TEA, buffered phosphate or bicarbonate
- Incompatible solvents flush the system well!
 - Avoid using DCM, Ethyl Acetate, DMSO, THF, MtBE and Chloroform with coated CSPs.
 - Okay with immobilized



Preemptive Retention or Selectivity Issue Solutions



Before connecting the column to the system: Flush all the HPLC unit with a compatible solvent – preferably 2-propanol. 1. Flush the entire unit with the column storage mobile phase. 2. Injector to column connection **Injection Loop** Pump to Detector **ALL** inlet solvent lines Injector (0010) **Injection Wash** connection Solvent -ee pf **Column to detector connection** Pump

Pressure



Preemptive Retention or Selectivity Issue Solutions







Regeneration or column flushing

- Repack the column with new CSP
- Chemical Attack

Adsorbed sample

Replace the column

©2022 Daicel Chiral Technologies. All rights reserved.



After Retention or Selectivity Changes Have Occurred











Services for Prep Columns

- Column Hardware Replacement
 - Open column, replace frit, ensure the CSP bed is not damaged
- Column Recondition
 - Take the column through a regeneration sequence remove and clean-up CSP, add in new CSP, repack, and QC

Less cost than purchasing a new column







www.chiraltech.com

©2022 Daicel Chiral Technologies. All rights reserved.





Table 1. Immobilized Primary Screening Solvents

<i><i>D</i></i> ^[']			Contact Us	How to Order	News & Events	Daicel Life Scien	ices
Chiral Technologies	SOLUTIONS	CHIRAL SELECTORS	CHIRAL TECHNI	QUES SERV	VICE & SUPPORT	ABOUT US	Q

Searchable Product Index

find the column you are looking for, please contact us via the link below.

Find the ideal column for your system using our searchable product index. Use the pull-down

Access Instruction Manuals by entering your Part Number and clicking the description link.

menus to narrow your selections from 1200+ Daicel chromatography columns. If you are unable to

Primary solvent mixtures	Alkane 0/2-PrOH	Alkane 0 /EtOH	Alkane 0 /MtBE/EtOH 0	Alkane 0 /THF 8	Alkane/DCM G /EtOH
Typical starting conditions	80:20	80:20	0:98:2	70:30	50:50:2
Advised optimisation range	99:1 to 50:50	99:1 to 50:50	80:20:0 to 0:40:60	95:5 to 0:100	85:15:0 to 0:80:20

• Alkane = n-Hexane, iso-Hexane or n-Heptane. Some small selectivity differences may sometimes be found.

In absence of alkane, methanol is more efficient than ethanol when combined with MBE.

In the case of no environmental restrictions, use of DCM is preferred to THF in terms of better enantioselectivity that the former may induce.
 For excessively retained samples, addition of ethanol up to 20% in pure DCM would be helpful.

D – Additives

For basic or acidic samples, it is necessary to incorporate an additive into the mobile phase in order to optimise the chiral separation.

• It has been found that certain amines, such as EDA and AE induce much better behaviour for certain basic compounds than the most commonly used DEA.

Basic Samples	Acidic Samples	
require	require	
Basic additives	Acidic additives	
Diethylamine (DEA) 2-Aminoethanol (AE) Ethylenediamine (EDA) Butyl amine (BA)	Trifluoroacetic acid (TFA) Acetic acid Formic acid	
< 0.5%	< 0.5%	
Typically 0.1%	Typically 0.1%	

☞ The addition of a low percentage of an alcohol (e.g. 2% EtOH or MeOH) in the mobile phase may be helpful to ensure the miscibility of EDA and AE with the low polarity mobile phases.

STRONGLY BASIC solvent additives or sample solutions <u>MUST BE AVOIDED</u>, because they are likely to damage the silica gel used in this column.

Column Care / Maintenance

- **D** The use of a guard cartridge is highly recommended for maximum column life.
- Samples should be dissolved in the mobile phase. The mobile phase and the sample solution should be filtered through a membrane filter of approximately 0.5µm porosity to ensure that there is no precipitate before using.

Column cleaning and regeneration procedures

Following extensive use of the column in multiple solvents there may be a change in column reproducibility. In order to ensure consistent performance, a regeneration method may be implemented to eliminate any change in chiral recognition due to the history of the column (mobile phases, additives...).

• Flush with ethanol at 0.5 ml/min for 30 min, followed by 100% THF at 0.5 ml/min for 2 hours.

Help?

• Flush with ethanol at 0.05 ml/min^(*) for 300 min.



PART NUMBER	
DAICEL STATIONARY PHASE	TECHNIQUES
SELECTORS	PARTICLE SIZE (µm)
I.D. (mm)	LENGTH (mm)

©2022 Daicel Chiral Technologies. All rights reserved.









Pressure >



In Summary – Major Causes of Column Failure



• Column or inlet frit fouling

- Can be avoided with good sample cleanup
- Filters or guard columns
- Bed compression void formation
 - Operate at low or reasonable flow rate
 - Avoiding mechanical or thermal "shock"
- Chemical attack or modification
 - Use appropriate pH range mobile phases
 - Avoid restricted solvents for coated phases







• <u>questions@cti.daicel.com</u> – General column chromatography or method questions

• <u>wumstead@cti.daicel.com</u> – Presentation, chromatography, or column questions



