



CHIRAL SELECTORS NAMES & STRUCTURES

STRUCTURES OF COMMERCIALY AVAILABLE DAICEL CHIRAL STATIONARY PHASES

AMYLOSE DERIVATIVES











CELLULOSE DERIVATIVES















IMMOBILIZED POLYSACCHARIDE SELECTORS

PRIMARY SCREENING SET





COATED POLYSACCHARIDE SELECTORS





PROTEIN-BASED SELECTORS

CHIRALPAK AGP ●●●●● α ₁ -acid glycoprotein	CHIRALPAK BSA ●●●●● bovine serum albumin	CHIRALPAK CBH ●●●●● cellobiohydrolase (stable enzyme)	CHIRALPAK DSA ●●●●● dog serum albumin
CHIRALPAK HSA ●●●●● human serum albumin	CHIRALPAK MSA ●●●●● mouse serum albumin	CHIRALPAK PSA ●●●●● pig serum albumin	CHIRALPAK RSA ●●●●● rat serum albumin





SELECTORS FOR SPECIAL APPLICATIONS



ZWITTERION EXCHANGER





ANION EXCHANGER



LIGAND EXCHANGER



CROWN ETHER AND IMMOBILIZED CROWN ETHER



LEGEND

- <2 μm
- 3 μm
- 5 μm
- 10 μm
- 20 μm

USP CODE DESIGNATIONS

L32 CHIRALPAK WH	L40 CHIRALCEL OD CHIRALCEL OD-H	L41 CHIRALPAK AGP	L51 CHIRALPAK AD CHIRALPAK AD-H CHIRALPAK AD-3	L66 CROWNPAK CR(+)
L70 CHIRALCEL OC-H	L79 CHIRALPAK HSA	L80 CHIRALCEL OJ CHIRALCEL OJ-H CHIRALCEL OJ-3	L90 CHIRALPAK AS CHIRALPAK AS-H CHIRALPAK AS-3	L93 CHIRALCEL OD-RH CHIRALCEL OD-3R
L99 CHIRALPAK IA CHIRALPAK IA-3	L107 CHIRALCEL OJ-RH	L108 CHIRALPAK CBH	L117 CROWNPAK CR(-)	

CHIRAL TECHNOLOGIES / DAICEL GROUP – LOCATIONS

NORTH/LATIN AMERICA

Chiral Technologies, Inc.
800 North Five Points Road
West Chester, PA 19380
Tel: +1-610-594-2100
Fax: +1-610-594-2325
Email: chiral@chiraltech.com
www.chiraltech.com

CHINA

Daicel Chiral Technologies (China) Co., Ltd
No. 32, HeXiang Road
Waigaoqiao FreeTradeZone, Post code: 200131
Shanghai, China
Tel: +86 21 504 60086
Fax: +86 21 504 62321
Email: chiral@cdtc.daicel.com
www.daicelchiraltech.cn

EUROPE

Chiral Technologies Europe SAS
Parc d'Innovation -
160, Bd Gonthier d'Andernach - CS 80140
67404 Illkirch Cedex FRANCE
Tel: +33 (0) 3 88 79 52 00
Fax: +33 (0) 3 88 66 71 66
Email: cte@chiral.fr
www.chiraltech.com

ASIA

Daicel Corporation
CPI Company
JR Shinagawa East Building
2-18-1, Konan, Minato-ku
Tokyo 108-8230 JAPAN
Tel: +81 (0) 3 6711 8221
Fax: +81 (0) 3 6711 8228
Email: chiral@jp.daicel.com
www.daicelchiral.com

INDIA

Daicel Chiral Technologies (India) Pvt. Ltd
Lab No: A4, Phase III, IKP Knowledge Park
Genome Valley, Turkapally, Shameerpet
Ranga Reddy District
Hyderabad - 500 078, Telangana INDIA
Tel: + 91 40 2338 3700 or 91 40 2348 0103
Fax: +91 40 2348 0104
Email: chiral@chiral.daicel.com
www.chiraltech.com



WWW.CHIRALTECH.COM

CHIRALCEL, CHIRALPAK and CROWNPAK are registered trademarks of DAICEL CORPORATION

