

Chiral Technologies

USP List of HPLC Columns



USP Code Column Name	USP Specification Analysis (Pending/Official)	Drug Chemical Name	Drug Brand Name
L32 CHIRALPAK® WH	Limit of Organotin	Fluorodopa	FDOPA
L40 CHIRALCEL® OD CHIRALCEL OD-H CHIRALCEL OD-3	Limit of efavirenz R-enantiomer Limit of voriconazole-related compound B Limit of atomoxetine S- enantiomer Limit of valsartan-related compound A Limit of irinotecan hydrochloride enantiom	Efavirenz Voriconazole Atomoxetine Hydrochloride Valsartan Irinotecan Hydrochloride	Sustiva VFEND Strattera Diovan Camptosar
L41 CHIRALPAK AGP	Limit of ropivacaine-related compound B Limit of galantamine 4R, 8R-stereoisomer Limit of montelukast S-enantiomer Limit of tenofovir disoproxil-related compo	Ropivacaine Galantamine Hydrobromide Montelukast Sodium Tenofovir Disoproxil Fumarate	Naropin Razadyne Singulair Viread
L51 CHIRALPAK AD CHIRALPAK AD-H CHIRALPAK AD-3	Limit of abacavir enantiomer Limit of levetiracetam R-enantiomer Limit of tolterodine S-enantiomer Limit of atorvastatin 3S,5S-steriomer Measurement of diastereoisomeric ratio Limit of paroxetine-related compound C Limit of brinzolamide-related compound A Enantiomeric Purity	Abacavir Sulfate Levetiracetam Tolterodine Tartrate Atorvastatin Calcium Fulvestrant Paroxetine Hydrochloride Brinzolamide Darifenacin Hydrobromide	Ziagen Keppra Detrol Lipitor Faslodex Paxil Azopt Enablex
L66 CROWNPAK® CR(+)	Limit of valacyclovir D-isomer	Valacyclovir Hydrochloride	Valtrex
L70 CHIRALCEL OC-H	Limit of oxaliplatin-related compound D	Oxaliplatin	Eloxatin
L75 CHIRALPAK BSA	Enantiomeric Purity	Methotrexate	Otrexup
L79 CHIRALPAK HSA	Enantiomeric Purity	Calcium L-5- Methyltetrahydrofolate	Calcium Supplements
L80 CHIRALCEL OJ CHIRALCEL OJ-H CHIRALCEL OJ-3	Limit of Exemestane To be specified	Exemestane Clopidogral Bisulfate	Aromasin Plavix

Chiral Technologies

USP List of HPLC Columns



USP Code Column Name	USP Specification Analysis (Pending/Official)	Drug Chemical Name	Drug Brand Name
L90 CHIRALPAK AS CHIRALPAK AS-H CHIRALPAK AS-3	Enantiomeric Purity	Frovatriptan Succinate	Frova
L93 CHIRALCEL OD-RH CHIRALCEL OD-3R	Organic Impurities	Ezetimibe	Zetia
L99 CHIRALPAK IA CHIRALPAK IA-3	Enantiomeric Purity	Darifenacin Hydrobromide	Enablex
L107 CHIRALCEL OJ-RH CHIRALCEL OJ-3R	To Be Specified	Rosuvastatin Calcium	Crestor
L108 CHIRALPAK CBH	Enantiomeric Purity	Epinephrine Bitartrat	EpiPen, Adrenalin
L117 CROWNSPAK CR(-)	Identification and Related compound	Lysine	Lysine
L119 CHIRALPAK IC CHIRALPAK IC-3	Organic Impurities	Ezetimibe and Simvastatin	Vytorin

LOCATIONS:

North/Latin America

Chiral Technologies, Inc.
800 North Five Points Road
West Chester, PA 19380
800 6 CHIRAL
Tel: 610-594-2100
Fax: 610-594-2325
chiral@cti.daicel.com
www.chiraltech.com

Europe

Chiral Technologies Europe
Parc d'Innovation
Bd Gonthier d'Andernach
67400 Illkirch Cedex, France
Tel: +33-388-795-200
Fax: +33-388-667-166
cte@cte.daicel.com
www.chiraltech.com

India

Daicel Chiral Technologies (India) Pvt. Ltd.
Survey No.542/2, IKP Knowledge Park,
Turkapally, Shamirpet Mandal,
Medchal-Malkajgiri District,
Hyderabad-500101. Telangana, India
Tel : +91 84 1866 0700 or 91 84 1866 0703
Fax : +91 84 1866 0730
Email: chiral@chiral.daicel.com
www.chiraltech.com